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Table 5.1.4.-1. – Acceptance criteria for microbiological quality of non-sterile dosage forms

Route of administration	TAMC (CFU/g or CFU/mL)	TYMC (CFU/g or CFU/mL)	Specified micro-organisms
Non-aqueous preparations for oral use	10 ³	10 ²	Absence of <i>Escherichia coli</i> (1 g or 1 mL)
Aqueous preparations for oral use	10 ²	10 ¹	Absence of <i>Escherichia coli</i> (1 g or 1 mL)
Rectal use	10 ³	10 ²	-
Oromucosal use Gingival use Cutaneous use Nasal use Auricular use	10 ²	10 ¹	Absence of <i>Staphylococcus aureus</i> (1 g or 1 mL) Absence of <i>Pseudomonas aeruginosa</i> (1 g or 1 mL)
Vaginal use	10 ²	10 ¹	Absence of <i>Pseudomonas aeruginosa</i> (1 g or 1 mL) Absence of <i>Staphylococcus aureus</i> (1 g or 1 mL) Absence of <i>Candida albicans</i> (1 g or 1 mL)
Transdermal patches (limits for one patch including adhesive layer and backing)	10 ²	10 ¹	Absence of <i>Staphylococcus aureus</i> (1 patch) Absence of <i>Pseudomonas aeruginosa</i> (1 patch)
Inhalation use (special requirements apply to liquid preparations for nebulisation)	10 ²	10 ¹	Absence of <i>Staphylococcus aureus</i> (1 g or 1 mL) Absence of <i>Pseudomonas aeruginosa</i> (1 g or 1 mL) Absence of bile-tolerant gram-negative bacteria (1 g or 1 mL)
◆ Special Ph. Eur. provision for oral dosage forms containing raw materials of natural (animal, vegetal or mineral) origin for which antimicrobial pretreatment is not feasible and for which the competent authority accepts TAMC of the raw material exceeding 10 ³ CFU/g or CFU/mL.	10 ⁴	10 ²	Not more than 10 ² CFU of bile-tolerant gram-negative bacteria (1 g or 1 mL) Absence of <i>Salmonella</i> (10 g or 10 mL) Absence of <i>Escherichia coli</i> (1 g or 1 mL) Absence of <i>Staphylococcus aureus</i> (1 g or 1 mL)◆
◆ Special Ph. Eur. provision for premixes for medicated feeding stuffs for veterinary use using excipients of plant origin for which antimicrobial treatment is not feasible.	10 ⁵	10 ⁴	Not more than 10 ⁴ CFU of bile-tolerant gram-negative bacteria (1 g or 1 mL) Absence of <i>Escherichia coli</i> (1 g or 1 mL) Absence of <i>Salmonella</i> (25 g or 25 mL)◆

Microbiological monitoring. During production and subsequent storage, appropriate measures are taken to ensure that the microbial count is adequately controlled and monitored. Appropriate alert and action levels are set so as to detect adverse trends. Under normal conditions, an appropriate action level is a microbial count of 10 CFU per 100 mL when determined by filtration through a membrane with a nominal pore size not greater than 0.45 µm, using R2A agar, using at least 200 mL of water for injections in bulk and incubating at 30-35 °C for not less than 5 days. For aseptic processing, stricter alert levels may need to be applied.

EP Monograph on water for injection

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*"Unlike non-sterile dosage forms, there are **no recommended bioburden levels provided in regulatory guidelines or compendia for the [downstream]** protein purification processes of biologic or other biopharmaceutical products, therefore, manufacturers are responsible for setting bioburden control levels for biologic production processes."*

"The BPOG Bioburden Working Group conducted a member survey of bioburden action levels and found that the majority of biologic processes action levels were set between 1-10 CFU/mL."

Bain, D. "Microbial Monitoring For Biological Drug Substance Manufacturing: An Industry Perspective"
BioPhorum Operations Group. 2015

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

- 目標
- EP 2.6.12 and USP <61>: 'The tests described will allow quantitative enumeration of mesophilic bacteria and fungi that may grow under aerobic conditions. The tests are designed to determine whether a substance or preparation complies with an established specification for microbiological quality.'
- “所描述的測試將允許定量計數可能在有氧條件下生長的嗜溫細菌和真菌。這些測試旨在確定物質或製劑是否符合既定的微生物品質規範。”

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

- 目標:
- EP 2.6.12 and USP <61>:
- 'Test each batch of ready-prepared medium and each batch of medium prepared either from dehydrated medium or from the ingredient described'
- '測試每批現成培養基和每批由脫水培養基或所述成分製備的培養基'

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	Media	Test strain	Culture condition
EP8.3	R2A agar	Pseudomonas aeruginosa ATCC 2027	30 - 35°C Within 3day
		Bacillus subtilis ATCC 6633	
		Preculture by SCD、dilute by NaCl peptone buffer、inoculate <100CFU, use this strain solution within 2hours, use them at 2 ~ 8°C within 24hours (for test sample also)	
USP43 (1231)	Plate Count Agar Or TGYA medium as high nutrient R2A as low nutrient	No concrete information about standard strain and culture condition=>changed.. <i>Pseudomonas aeruginosa ATCC 9027 and Bacillus subtilis ATCC 6633.</i> <i>Additional organisms should be used to represent those that are considered objectionable and/or typically isolated from the water system (house isolates).</i>	30 - 35°C Within 3day <i>(Or low temp at least 4days)</i>
CHP 2015	R2A	Pseudomonas aeruginosa CMCC(B)10 104	30 - 35°C at least 5 days
		Bacillus subtilis CMCC(B) 63 501	
		Preculture by TSB or TSA at 30-35°C 18-24h, inoculate ≤100CFU	

- “The type of filter material is chosen such that the bacteria retaining efficiency is not affected by the components of the sample to be investigated.”
- “選擇過濾材料的類型，以確保細菌截留效率不受待研究樣品成分的影響。”

EP 2.6.12 Microbial Examination of non-sterile products

- Select the material of the membrane depending on the test product (e.g. antibiotics: use a low absorption membrane)
- “Adherence of residual antimicrobial agents to the filter membrane may cause growth inhibition. Filtration through a low-binding filter material, such as polyvinylidene difluoride, helps to minimize this growth inhibition..”
- 根據測試產品選擇膜的材料（例如抗生素：使用低吸收膜）
- “殘留的抗菌劑粘附在濾膜上可能會導致生長抑制。通過低結合力過濾材料（如聚偏二氟乙烯）進行過濾，有助於最大限度地減少這種生長抑制。

USP <1227> validation of microbial recovery from pharmacopeial articles

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



USP <61> & USP <62> Microbiological examination of non sterile products

EP 2.6.12 Microbial examination of non sterile products

USP <71> Sterility testing

USP <1231> Water for Pharmaceutical purposes

Choose your membrane wisely 明智地選擇您的濾膜

Membrane Filtration: Retention & Recovery

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Retention 保留


- Ability to capture all microorganisms on membrane surface
能夠捕捉膜表面的所有微生物

Recovery 回收

- Ability to grow captured microorganisms
- Proper nutrients, temperature, and time
- Filter matrix must allow free flow of nutrients
- Cell must be able to replicate
- 能够生长捕获的微生物
- 适当的营养物质、温度和时间过滤基质必须允许营养物质的自由流动
- 细胞必须能够复制

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Examined relationship of membrane pore size to growth characteristics and recovery on membranes used for enumeration tests
檢測膜孔徑與用於列舉試驗的膜的生長特性和恢復的關係

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- 0.2μm, 0.45μm 0.7μm, 0.8μm, 0.65μm, 1.2μm

Examined a variety of microorganisms
檢測了各種微生物

B. subtilis 13933

C. albicans 10231

E. aerogenes 11437

M. luteus 9341

Primary effluent

Br. diminuta 19146

C. sporogenes 11437


E. coli 25922

P. agglomerans (environmental)

Pore size considerations for bioburden testing method

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- 0.45 μm retained the microorganisms as well as 0.22 μm
- 0.45 μm consistently demonstrated acceptable recovery (>90% vs spread plate) of a broad range of microorganisms
- 0.22 μm did not always demonstrate acceptable recovery
- Colony morphology was similar to or slightly smaller for 0.45 μm and 0.22 μm vs spread plates but still acceptable
- Standard Methods currently recommend 0.45 μm membranes for microbiology testing

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EU GMP Annex 1 2022

9.30 **Action limits** for viable particle contamination are shown in Table 6

Table 6: Maximum action limits for viable particle contamination

Grade	Air sample CFU /m ³	Settle plates (diam. 90 mm) CFU /4 hours ^(a)	Contact plates (diam. 55mm), CFU / plate ^(b)	Glove print, Including 5 fingers on both hands CFU / glove
A	No growth ^(c)			
B	10	5	5	5
C	100	50	25	-
D	200	100	50	-

(a) Settle plates should be exposed for the duration of operations and changed as required after 4 hours. Exposure time should be based on recovery studies and should not allow desiccation of the media used.

(b) It should be noted that for Grade A, the expected result should be no growth.

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United States Pharmacopeia (USP) <1116> Microbiological evaluation of clean rooms and other controlled environments May 2012

Suggested initial Microbial Contamination

Recovery rates (%) in Aseptic environments (Targeted Levels)

Room Classification	Active Air Sample (%)	Settle Plate (9 cm) 4h Exposure (%)	Contact Plate or Swab (%)	Glove or Garment (%)
Isolator/Closed RABS (ISO 5 or better)	<0.1	<0.1	<0.1	<0.1
ISO 5	<1	<1	<1	<1
ISO 6	<3	<3	<3	<3
ISO 7	<5	<5	<5	<5
ISO 8	<10	<10	<10	<10

(a) All operators are aseptically gowned in these environments (with the exception of background environments for isolators). These recommendations do not apply to production areas for non-sterile products or other classified environments in which fully aseptic gowns are not donned.

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

5.5.2 Culture media

Culture media shall, if not otherwise indicated, be non-selective.

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International
Organization for
Standardization

- ISO 14698-1

5.5.3 Incubation

Total incubation periods of two to five days for bacteria and five to seven days for fungi are generally acceptable, especially when the number of VU is low. When anaerobic, thermophilic, micro-aerophilic, or nutritionally deficient or fastidious bacteria, and fungi are of concern, specific atmospheric conditions and incubation times could be necessary. Plates should be observed at appropriate intervals over the incubation period.

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The microbiological culture media used in environmental monitoring should be validated as capable of detecting fungi (i.e., yeasts and molds) as well as bacteria and incubated at appropriate conditions of time and temperature. Total aerobic bacterial count can be obtained by incubating at 30 to 35°C for 48 to 72 hours. Total combined yeast and mold count can generally be obtained by incubating at 20 to 25°C for 5 to 7 days.

Selection of Culture Conditions

Incubation time and incubation temperatures are set once the appropriate media have been selected. Typically, for general microbiological growth media such as SCDM, incubation temperatures in the ranges of approximately 20°–35° have been used with an incubation time of not less than 72 hours. Longer incubation times may be considered when contaminants are known to be slow growing. The temperature ranges given above are by no means absolute.

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ZHEJIANG FANGYUAN TEST GROUP CO., LTD.